

HEALTHY


- ▶ HEALTHY GINGIVA
- ▶ NO BLEEDING ON PROBING



NOT NEEDED



PS TIP IF CALCULUS



GINGIVITIS

- ▶ SHALLOW POCKET
- ▶ NO BONE LOSS
- ▶ BLEEDING ON PROBING




4 mm

NOT NEEDED



PS TIP IF CALCULUS



PERIODONTITIS

- ▶ DEEP POCKET
- ▶ BONE LOSS
- ▶ BLEEDING ON PROBING
- ▶ NO SUPPURATION




4 mm



9 mm



10 mm PS tip



SUPPURATING PERIODONTITIS

- ▶ BONE LOSS
- ▶ BLEEDING ON PROBING
- ▶ SUPPURATION



ONLY SUPRA



10 mm PS tip